


KN95 MASKS



PRODUCT INFORMATION

TYPE:	KN95 Surgical Mask
COMPONENT:	Non-woven cloth with an ear loop mask body, nasal bridge piece, and a mask belt
FUNCTION + APPLICATION:	Respirator mask worn on the face used to block the exhalation of the oral and nasal cavities or to eject pollutants in a medical environment.
CAUTION	(1) Do not use when expired (2) A one-time use only (3) Use with caution if allergic to non-woven material (4) Refer to instructions before use
STORAGE	Store at room temperature in a non-corrosive environment with no more than 80% humidity. Avoid high temperatures.

 www.atlanticpkg.com

 (800) 722-5841

